

New Design. Touch Screen. Multiple Analyte Reading & More.





Ready-to-use ELISA instrument programmed with Monobind assays and proprietary formulas designed, optimized and certified by Monobind for two-point analysis.

> Features

- Standalone ELISA reader with integrated computer and printer
- Easy and intuitive touch-screen operation with stylus
- Dynamic reading to 4.5 OD or above (from 3.0 maximum with Eldex 3.8)
- Carrier with (new) automated movement handles 3 strips with 1 to 12 wells/ strip
- Ability to program and read different analytes in strip 1 through 3 or to read all 3strips with the same analyte automatically
- Full 96 well-plate read in 3 to 4 shifts (2 minutes per load/ shift)
- Test storage increased from 100 to 120 tests
- Ten digit alphanumeric sample ID now available for patient tracking
- Password protection system enables open or closed instrument use
- Administrator tool provides summary of wells read per each stored test to assist distributors with laboratory-use monitoring
- USB port easily imports new software or exports readings, including graphs to a PC for archiving
- 110/220 power options, light weight, small footprint
- Optics include 4 filters (405, 450, 492 and 630)
- Quality instrument, US-built, CE mark and ISO standards, NRTL certified



Broad Spectrum of 30+ Analytes

NeoEldex is completely integrated with our immunoassays

- Programmed with Monobind's 50+ AccuBind® ELISA tests, parameters and data analytics for optimal results
- Exceptionally easy: Simply select test using touch screen then insert wells
- Exclusive 2PT Calibration enables labs to cut cost-per-test by using two calibrators (instead of six or seven) on same lot kits within 30-days of running full curve
- Run with Monobind multi-ligand controls to ensure total quality, recommended for 2PT calibration



Available in 45+ Analytes

Contact us today to learn more

www.monobind.com :: sales@monobind.com

